



BETTER BASICS
LABORBEDARF

1ST QUARTER 2026




Modules



Product catalogue

Holder modules

Modules for SmartRack® and SmartIntegrate



www.Better-Basics-Laborbedarf.de

Unlock Your Lab Potential

Table of Contents

Your Advantage in Everyday Laboratory Work.....	3
Modules I for SmartRack®- and SmartIntegrate-Products.....	6
The Configuration.....	6
Modules.....	7
Additional accessories.....	30
SmartStand.....	31
Customising.....	32
Item Overview.....	34
Technical Information.....	38

Your Advantage in Everyday Laboratory Work

Better Basics Laborbedarf – Your Advantage in Everyday Laboratory Work

Better Basics Laborbedarf optimises laboratories and their workspaces with innovative, award-winning organisational systems. We support our customers in creating the best possible working environment for their laboratory staff, helping to unlock the full potential of their laboratories. This approach ensures consistently high quality, efficiency, and process stability in the lab while providing a safe, ergonomic workplace with satisfied employees.

Enhancing Employee Satisfaction

Our SmartRack® and SmartIntegrate product families offer a comprehensive system to optimise laboratory workflows. With over 100 modular solutions, they not only ensure order but also noticeably improve working conditions. Enhanced ergonomics, reduced movement paths, and clearly structured workstations increase both efficiency and employee satisfaction. Safe and organised environments reduce stress and promote job satisfaction. Furthermore, they demonstrably contribute to improved workplace safety and occupational health through thoughtful placement and ease of use.


A System for Sustainable Process Improvement

Our product families are designed to standardise workflows and save time. The modular concept of SmartRack® and SmartIntegrate allows for flexible adaptations to any laboratory environment – from safety workbenches to large open-plan labs. These systems not only optimise space utilisation but also increase efficiency by an average of 30% and save up to 80% space. With a wide range of modules and accessories, they offer customised solutions for every need.

Efficient, Safe, and Satisfying

For decision-makers and managers, SmartRack® and SmartIntegrate ensure process stability and economical resource use. At the same time, these systems provide innovation managers and organisational developers with the ideal basis for forward-looking workplace designs. For laboratory staff, this means fewer errors thanks to well-structured work areas, a more relaxed working environment, and often the ability to finish earlier due to significantly improved workflows.

Discover how Better Basics Laborbedarf can transform your everyday laboratory work – for more efficient processes, satisfied employees, and excellent results.



The product overview of Better Basics

▶ SmartRack®, the mobile solution for: laboratory workbench, fume hood and safety cabinet



+ Customised solutions

▶ SmartIntegrate, the integrative solutions for

I. Laboratory workbench



II. Fume hood



III. Safety cabinet



IV. Integration into laboratory equipment



+ Customised solutions

▶ Modules

- Currently, more than 100 modules
- All modules in the Better Basics catalogue can be used universally in all SmartRack® and SmartIntegrate solutions
- Customised modules can be developed based on the specific customer needs, all customised modules can be used in all SmartRack® and SmartIntegrate solutions



+ Customised solutions

Our Modules - Versatility for Every System

The modules from Better Basics are the ideal building blocks for a well-organised laboratory. They can be seamlessly integrated into the SmartRack® system as well as the Smart-Frame, SmartRail, and SmartEdge solutions from the SmartIntegrate product line.

Whether pipette holders, storage boxes, or inserts for specialised applications, all modules are easy to switch and ensure precise adaptation to your workflows. Their innovative design and robustness contribute to the efficiency and ergonomics of every workstation.



Photo: A selection of innovative modules for the SmartRack® and SmartIntegrate product line

Modules | The configuration

<p>1. Selection of organisational solution:</p> <p>SmartIntegrate SmartRail, SmartEdge, SmartFrame</p>			
<p>SmartRack® SmartRack®, SmartRack® mini</p>			
<p>2. Selection of the modules: Holders and storage box and storage box inserts</p>			



Photo: Two SmartRack® and one SmartFrame MB 001 with a SmartStand in a Wesemann fume hood

Modules

Storage box



The module is the basis for storing and fixing various smaller laboratory equipment. The storage box can hold a maximum of three storage box inserts, which can be arranged and combined as required. Alternatively, the storage box can also be used without storage box inserts.

Item	Type	PU	Item no.	Price in EUR
Storage box	for 2 storage box inserts	1	01185	69,99
Storage box	for 3 storage box inserts	1	01142	79,99

Storage box insert – for small accessories



The module can be inserted as required into the storage box and enables easy storage and fixation of small accessories (e.g. pens, spatulas or tweezers).

Storage box and small accessories are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Storage box insert – for small accessories		1	01148	64,99

Storage box insert – for microreaction vessels



The module can be inserted as required into the storage box and enables easy storage and fixation of microreaction vessels.

Storage box and microreaction vessels are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Storage box insert – for microreaction vessels	32 x 0,5 ml	1	01104	64,99
Storage box insert – for microreaction vessels	20 x 2 ml	1	01130	64,99

Storage box insert – for sample containers and test tubes



The module can be inserted as required into the storage box and enables easy storage and holding of sample containers and test tubes. It is equipped with a spring mechanism that provides a secure hold for sample vials and test tubes of different diameters.

Storage box, sample containers and test tubes are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Storage box insert – for sample containers and test tubes	20 x Ø 3 - 10 mm	1	01161	64,99
Storage box insert – for sample containers and test tubes	15 x Ø 8 - 12 mm	1	01181	64,99
Storage box insert – for sample containers and test tubes	12 x Ø 10 - 20 mm	1	01160	64,99
Storage box insert – for sample containers and test tubes	8 x Ø 18 - 22 mm	1	01182	64,99
Storage box insert – for sample containers and test tubes	6 x Ø 20 - 30 mm	1	01159	64,99

Storage box insert – for centrifuge tubes



The module can be inserted as required into the storage box and enables easy storage and fixation of centrifuge tubes.

Storage box and centrifuge tubes are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Storage box insert – for centrifuge tubes	12 x 15 ml	1	01128	64,99
Storage box insert – for centrifuge tubes	6 x 15 ml, 4 x 50 ml	1	01102	64,99
Storage box insert – for centrifuge tubes	6 x 50 ml	1	01129	64,99

Adapter for LAUDA lid holder



The adapter is used to hold a LAUDA cover holder for storing thermostat covers.

Item	Type	PU	Item no.	Price in EUR
Adapter for LAUDA lid holder		1	01223	49,99

Holder for 2/3 pipettes



The module allows easy attachment of up to three pipettes with hook. The pipette tips are not affected and therefore not contaminated.

Pipettes are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for 2 pipettes		1	01180	64,99
Holder for 3 pipettes		1	01145	69,99

Holder for 3 Eppendorf pipette holders



The module enables the safe attachment of a pipetting aid.

Pipettes are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for 3 Eppendorf pipette holders		1	01204	79,99

Holder for barcode scanner



The module enables easy attachment of a barcode scanner.

Barcode scanner is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for barcode scanner	suitable for scanners such as Newland NLS-HR22 and Aplic 304674	1	01163	74,99

Holder for cup



The module enables the safe storage of a 500 ml AKU® drinking cup, which can be used as storage for consumables such as plastic transfer- or disposable pipettes.

Cup is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for cup	for AKU® PP-Drinking cup, 500 ml	1	01213	59,99
AKU® PP-Drinking cup, 500 ml		1	09006	5,99

Holder for display



The module enables secure mounting of screens with a VESA 75-100 interface. The integrated ball joint (only for Item no. 01221-01222) allows flexible adjustment in all directions, ensuring optimal positioning while keeping the screen firmly in place.

Display is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for display	for VESA 75, ball joint	1	01221	229,99
Holder for display	for VESA 100, ball joint	1	01222	229,99
Holder for display	for VESA 75, fixed	1	01230	169,99
Holder for display	for VESA 100, fixed	1	01231	169,99

Holder for chromatography column



The module allows effortless attachment of chromatography columns. The bracket adjusts to the column circumference via a snap-in mechanism to provide a secure fit. If the length of the suspended equipment exceeds 15 cm, we recommend the use of additional holder modules, stacked in the required number, to ensure the stability of the set-up. Both the outside of the chromatography column and the silicone covers of the chromatography column holders must be clean and free of dust and grease.

Chromatography column is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for chromatography column	Ø 5 - 10 mm	1	01177	69,99
Holder for chromatography column	Ø 10 - 15 mm	1	01176	69,99
Holder for chromatography column	Ø 15 - 25 mm	1	01165	69,99
Holder for chromatography column	Ø 25 - 35 mm	1	01107	69,99
Holder for chromatography column	Ø 35 - 46 mm	1	01139	69,99
Holder for chromatography column	Ø 46 - 56 mm	1	01179	79,99

Holder for various KNAUER accessories



The module enables the safe attachment of various KNAUER accessories.

KNAUER accessories are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for various KNAUER accessories		1	01210	59,99

Holder for Erlenmeyer flasks



The module allows easy attachment of conical flasks with a beaded rim. The curved front means the conical flask is securely fixed in the module and is protected against falling out.

Erlenmeyer flasks are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for Erlenmeyer flasks	Erlenmeyer flasks wide neck 25, 100 ml, Erlenmeyer flasks narrow neck 250, 500 ml	1	01108	69,99
Holder for Erlenmeyer flasks	Erlenmeyer flasks wide neck 250 ml, 500 ml, 1000 ml, Erlenmeyer flasks narrow neck 2000 ml, 3000 ml, 5000 ml	1	01140	69,99

Holder for glassware with standard taper



The module allows easy attachment of glassware with standard ground joint ST14. The specially shaped recess prevents the glassware from falling out. Please note the information on standard ground joint tolerances of various standard taper (ST) on page 38.

Glassware is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for glassware with standard taper	ST 14	1	01111	49,99
Holder for glassware with standard taper	ST 29	1	01114	49,99

Holder for glassware with standard taper plus extension



The module corresponds in design and function to the module holder for glass devices with standard taper. However, the extended version is particularly suitable for easily hanging particularly large round-bottomed flasks or separating funnels, for example. Please note the information on standard ground joint tolerances of various standard taper (ST) on page 38.

Glassware is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for glassware with standard taper plus extension	ST 14	1	01113	59,99
Holder for glassware with standard taper plus extension	ST 29	1	01116	59,99

Holder for glassware with standard taper plus closure



The module corresponds to the module holder for glassware with standard taper, but is also equipped with a rotatable locking ring that protects glassware with standard taper against falling out. Please note the information on standard ground joint tolerances of various standard taper (ST) on page 38.

Glassware is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for glassware with standard taper plus closure	ST 14	1	01112	64,99
Holder for glassware with standard taper plus closure	ST 29	1	01115	64,99

Holder for glove / tissue boxes



The module allows easy attachment of glove or tissue boxes. The width of the bracket can be adjusted as required and adapted to the size of different boxes. Being able to pull out these consumables towards you makes it easy to remove them even when the box is low.

Tissue box is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for glove / tissue boxes		1	01143	89,99

Holder for Inlabtex Serial Dilution Bag box



The holder enables safe storage of an Inlabtex Serial Dilution Bag box.

Item	Type	PU	Item no.	Price in EUR
Holder for Inlabtex Serial Dilution Bag box		1	01216	79,99

Holder for cannula disposal container



The module enables the simple insertion of a cannula discharge box.

Cannula disposal container is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for cannula disposal container	for Vernacare sharpsafe® 1l	1	01174	49,99

Holder for KNAUER syringe filter box



The module enables the safe attachment of a KNAUER syringe filter box.

KNAUER syringe filter box is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for KNAUER syringe filter box		1	01212	79,99

Holder for KNAUER Valve Drive



The module enables the safe attachment of a KNAUER Valve Drive.

KNAUER Valve Drive is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for KNAUER Valve Drive		1	01217	79,99

Holder for kitchen roll



The module enables easy hanging of kitchen rolls. The roll can rotate in the module, ensuring single or multiple sheets can be unrolled easily.

Kitchen roll is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for kitchen roll		1	01162	89,99

Holder for cuvettes



The module makes it easy to store and fix cuvettes. Due to the special shape of the recesses, the side surfaces of the cuvette do not come into direct contact with the holder. The dimensions given refer to the layer thickness of the cuvette.

Cuvettes are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for cuvettes	10 x (10 mm x 10 mm)	1	01183	64,99

Holder for laboratory screw top bottles



The module enables easy attachment of laboratory screw top bottles. The module is equipped with a rotatable locking ring, which provides additional protection for the laboratory screw top bottles against falling out.

Laboratory screw top bottles are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for laboratory screw top bottles	10 - 25 ml (GL25)	1	01120	64,99
Holder for laboratory screw top bottles	50 ml (GL32)	1	01121	64,99
Holder for laboratory screw top bottles	100 - 1000 ml (GL45)	1	01122	64,99

Holder for microliter syringes



The module allows three microliter syringes to be safely attached. The snap-in mechanism provides a secure hold for the attached microliter syringes. The shape of the module shields the needles of the microliter syringes and thus protects the user against accidental contact with the tip of the needle.

Microliter syringes are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for microliter syringes		1	01155	89,99

Holder for bin liner




The module allows easy attachment of bin liners for small waste. Contamination can be prevented by disposing of waste directly where you are working.

The waste bag holder can optionally be fitted with a lid. This lid can be found in the "Accessories" section.

Bin liner is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for bin liner	without cover	1	01202	59,99

Holder for NETGEAR® ProSAFE GS108



The module enables the safe attachment of a NETGEAR® ProSAFE GS108.

NETGEAR® ProSAFE GS108 is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for NETGEAR® ProSAFE GS108		1	01211	69,99

Holder for notebook



The module allows a notebook to be placed on top. The integrated ball joint allows flexible adjustment in all directions, ensuring optimum positioning and a secure hold for the notebook.

Notebook is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for notebook	for 15 inch	1	01224	289,99

Holder for pipettes



The module enables easy attachment of pipettes. The pipette tips are not loaded and therefore not contaminated.

Pipettes are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for pipettes	for Eppendorf Multipipette® E3 / E3x / stram / Xstream	1	01189	59,99
Holder for pipettes	for Eppendorf Reference®	1	01164	49,99
Holder for pipettes	for Eppendorf Reference® 2	1	01168	49,99
Holder for pipettes	for Gilson® PIPETMAN® Classic / G	1	01166	49,99
Holder for pipettes	for INTEGRA VIAFLOW / INTEGRA VOYAGER	1	01184	49,99

Holder for pipette tip boxes



The module enables easy attachment of pipette tip boxes from different manufacturers and of different sizes. Due to a magnet integrated in the module and a magnetic stainless steel round inserted by the user into the pipette tip box, the secure hold of the pipette tip box is ensured. One stainless steel round is already included in the scope of delivery, additional stainless steel rounds can be ordered for working with multiple pipette tip boxes.

Pipette tip boxes are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for pipette tip boxes		1	01147	54,99
Stainless steel round for insertion into pipette tip boxes		1	01167	2,99

Holder for pipetting aid



The module enables the safe attachment of a pipetting aid.

Pipetting aid is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for pipetting aid	for Eppendorf Easypet® 3	1	01203	59,99
Holder for pipetting aid	for HIRSCHMANN® pipetus®	1	01188	59,99
Holder for pipetting aid	for INTEGRA-Pipetboy 2	1	01154	59,99
Holder for pipetting aid	for INTEGRA-Pipetboy pro	1	01156	59,99
Holder for pipetting aid	for BRAND accu-jet® pro	1	01178	59,99
Holder for pipetting aid	for BRAND accu-jet® S	1	01157	59,99
Holder for pipetting aid	for Thermo Scientific Matrix	1	01186	59,99
Holder for pipetting aid	for TPP Turbo-Fix	1	01207	59,99
Holder for pipetting aid	for VWR® Powerpette Pro	1	01199	59,99

Holder for separating funnel / shaking funnel



The module allows easy attachment of separating funnels / shaking funnels.

Note:

The diameter of the separating funnel / shaking funnel is measured at its widest point

Separating funnel is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for separating funnel / shaking funnel	Ø 45 - 85 mm	1	01152	59,99
Holder for separating funnel / shaking funnel	Ø 85 - 105 mm	1	01149	59,99
Holder for separating funnel / shaking funnel	Ø > 105 mm	1	01153	59,99

Holder for sensor systems



The module enables the safe attachment of a sensor system.

Sensor systems are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for sensor systems	for essentim sphere	1	01208	69,99
Holder for sensor systems	for essentim spot	1	01209	69,99
Holder for sensor systems	for essentim spider	1	01205	79,99

Holder for open fronted storage bins



The module enables the safe attachment of open fronted storage bins or stacking boxes with a mounting edge, which typically serves as a holding device for wall rails and shelving systems.

Open fronted storage bins are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for open fronted storage bins	suitable for storage bins such as SCHÄFER LF110, LF210, LF211, LF220	1	01158	69,99

Holder for wash bottles



The module enables easy attachment of wash bottles. The module can also be used to store screw top bottles of similar sizes.

Wash bottles are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for wash bottles	for 500 ml (up to \varnothing 75 mm)	1	01146	59,99
Holder for wash bottles	for 1000 ml (up to \varnothing 95 mm)	1	01173	69,99

Holder for tablet



The module allows easy attachment of tablets in portrait and landscape format. The tilted angle ensures good visibility and good usability of the display. The module can be adjusted to different tablet sizes via a locking mechanism, which can be adjusted in increments.

Tablet is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for tablet		1	01124	99,99

Holder for well plates



Well plates can be inserted into the module and also stacked. The horizontal orientation ensures optimal storage, e.g. of cell cultures.

Well plates are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for well plates		1	01125	54,99

Holder for ZEBRA printer ZD621



The module enables the safe storage of a ZEBRA printer ZD621.

Item	Type	PU	Item no.	Price in EUR
Holder for ZEBRA printer ZD621		1	01201	179,99

Holder for ZEBRA hand-held scanner TC22 standing at an angle



The module enables the safe storage of a ZEBRA hand-held scanner TC22 standing at an angle.

Scanner is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for ZEBRA hand-held scanner TC22 standing at an angle		1	01206	79,99

Holder for draining



The module allows smaller laboratory equipment to be placed on it to drain after rinsing.

Glassware is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Holder for draining	6x drainage prong, I-shaped	1	01225	59,99
Holder for draining	6x drainage prong, Y-shaped	1	01226	69,99
Holder for draining	3x drainage prong, I-shaped and 3x drainage prong, Y-shaped	1	01227	64,99


Shelf



The shelf can be used to store laboratory equipment, glassware and consumables up to a total weight of 3 kg, providing additional storage space for those items that cannot be stored in the SmartRack® support modules.

Item	Type	PU	Item no.	Price in EUR
Shelf	Width 299 mm, Depth 185 mm	1	01193	119,99

Dispenser for adhesive tapes



The tape or labeling tape dispenser is designed to hold up to 6 rolls with a core hole diameter of approximately 25 mm (1 inch) or approximately 76 mm (3 inches) with a total width of up to 160 mm. It is made of powder-coated steel, polyamide 12 and stainless steel. The serrated cutting device for easy cutting of the adhesive tape or labeling tape facilitates one-handed handling in the laboratory and a clean cut. The dispenser for adhesive tape or labeling tape is easy to clean. Two cores are included to accommodate rolls with a core hole diameter of approx. 25 mm (1 inch) or approx. 76 mm (3 inches).

Adhesive tapes are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Dispenser for adhesive tapes		1	01172	149,99
Additional core for adhesive tape dispenser	approx. 25 mm (1 inch) / approx. 51 mm (2 inch)	1	01198	14,99
Additional core for adhesive tape dispenser	approx. 25 mm (1 inch) / approx. 76 mm (3 inch)	1	01175	24,99

Dispenser for Parafilm® rolls



The dispenser for Parafilm® is designed to hold two standard rolls (50 mm) or one double roll (100 mm). It is made of aluminum and polyamide 12 and offers the possibility of clean and safe storage as well as easy removal of Parafilm® sealing film. Item no. 01214 is characterised by an integrated cutting device that enables the Parafilm® to be torn off.

Parafilm® rolls are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Dispenser for Parafilm® rolls	without tear-off blade	1	01171	149,99
Dispenser for Parafilm® rolls	with tear-off blade	1	01214	199,99

Dispenser for glass pasteur pipettes



The dispenser for glass Pasteur pipettes is made out of aluminum and polyamide 12 and polypropylene, has a front level indicator and holds up to 300 unpacked Pasteur pipettes of a length of 150 mm or 230 mm. The dispenser is easy to clean and offers the possibility of hygienic storage and individual removal of the pasteur pipettes from the dispenser. Filling is performed via a hinged lid.

Pasteur pipettes are not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Dispenser for glass pasteur pipettes	150 mm	1	01169	149,99
Dispenser for glass pasteur pipettes	230 mm	1	01170	169,99

Stand rod adapter



The stand rod adapter is hooked into the slotted rail like a module and offers an additional possibility to secure the SmartRack®. The supplied stand rod can be inserted into the stand rod adapter from all three directions in space and then fixed to an existing stand rod, or other permanently installed rods in the laboratory. This immobilises the SmartRack® and fixes it in that location.

Stand rod is not included in the scope of delivery.

Item	Type	PU	Item no.	Price in EUR
Stand rod adapter	ø 10 mm	1	01218	79,99
Stand rod adapter	ø 12 mm	1	01126	79,99

Additional Accessories

Lid for storage box



The lid for the storage box is an optional extension to the storage box. Please note that upgrading the storage box with a lid is only possible for the storage box from order year 2024 or later.

Item	Type	PU	Item no.	Price in EUR
Lid for storage box	for 2 storage box inserts	1	01229	34,99
Lid for storage box	for 3 storage box inserts	1	01228	49,99

Lid for bin liner holder



The lid for the bin liner holder (01202) is an optional extension to the bin liner holder. Please note that it is only possible to add a lid to the bin liner holder for bin liner holders ordered from 2024 onwards.

Item	Type	PU	Item no.	Price in EUR
Lid for bin liner holder		1	01192	19,99

SmartStand

SmartStand for well plates



The SmartStand for Well Plates can be rotated 360°, allowing the well plates to be filled from any side and the alignment to be adjusted quickly. The incrementally adjustable tilt of up to 40° makes the wells much easier to see from the user's perspective, which supports precise and ergonomic work when pipetting.



*SmartStand with well plate,
Well plate is not included in the scope of delivery.*



SmartStand without well plate.

Item	Type	PU	Item no.	Price in EUR
SmartStand for well plates		1	09001	129,99

Customising

Better Basics Laborbedarf offers an unparalleled approach to laboratory workstation organisation with its customisable and modular solutions. Below is an overview of how our customisation services span both the SmartRack® and SmartIntegrate product lines:

1. Customisation of the SmartRack®

For the SmartRack®

Our SmartRack® system is fully customisable to meet the specific requirements of your workspace. The dimensions of the SmartRack® can be individually adapted to the size of your equipment and workflows, ensuring seamless integration into your laboratory workstation. At the same time, the modular compatibility and reliability of your tailor-made SmartRack® are guaranteed through appropriately dimensioned slot rails. From compact configurations to larger frames, we enhance the organisation and efficiency of your laboratory.

2. Customisation of SmartIntegrate Products

a) Bespoke SmartFrame for Fume Hoods

The SmartFrame optimises the organisation within fume hoods, offering a structured and safe arrangement for laboratory glassware, tools, and equipment. Custom-built SmartFrames provide a perfect fit for your specific workspace.

b) Tailored SmartRail for Lab Benches

SmartRails bring adaptability and precision to your laboratory benches. Designed for both high-throughput workflows and compact setups, they ensure reduced clutter and enhanced efficiency. Each SmartRail is tailored to your preferred length and mounting specifications.

c) Custom Solutions for Safety Cabinets

For safety cabinets, our custom solutions include integrated SmartRails, SmartFrames and SmartEdges designed to maximise space utilisation while supporting your workflows. Each component is built to ensure easy cleaning and adherence to safety regulations.

3. Custom Modules

Tailored for any Laboratory Equipment

Using advanced 3D printing technology, we create modules and holders in any shape or size to fit your laboratory equipment. From pipette holders to mounts for specialised glassware, our bespoke designs cater to your unique requirements. Our rapid prototyping process ensures swift delivery of tailored solutions.

Why Choose Better Basics for Customisation?

- 1. Precision Craftsmanship:** Each component is tailored to your exact specifications, guaranteeing optimal performance and seamless integration.
- 2. Fast Delivery:** Our in-house design and manufacturing teams ensure swift turnaround from concept to delivery.
- 3. Transparent Pricing:** Costs for customisation are agreed upon in advance, providing clarity and control.
- 4. Scalability:** Whether expanding an existing laboratory or building a new one, our products adapt to meet your growing needs.

Transform Your Laboratory Workspace

With Better Basics Laborbedarf, your laboratory can benefit from customised SmartRack® or SmartIntegrate products.

Visit customising.better-basics-laborbedarf.de or contact our team today to unlock the full potential of your laboratory space.

This concise yet comprehensive approach ensures your laboratory benefits from enhanced efficiency, organisation, and safety—all tailored to your specific needs.

Examples for realised customisations

Special frame sizes

SmartRack® frame - customisable in width as well as height



SmartRack® frame - extra wide

For a special requirement we manufactured an extra wide frame of the SmartRack®. Due to the width, this SmartRack® makes optimal use of the working space in our customer's laboratory.

Special module sizes

Holder for multiple glassware



Double holder for glassware with ST10

For a project, a customer needed a double holder for glassware ST10. The module allows glassware with standard ground joint ST10 to be easily hung next to each other. The specially shaped recess prevents the glassware from falling out.

Of course, adaptations of this kind are also conceivable for other modules, e.g. storage boxes in individual sizes.

Special laboratory equipment

Dispenser for cylindrical Petri dish multi-packs



For a globally operating biotechnology company, a dispenser for cylindrical Petri dish multi-packs was developed.

The dispenser not only enables easy removal of the packs but also provides a quick overview of the remaining stock.

Item Overview

Modules

Item	Type	Item no.	PU	Price in EUR	Page
Storage box	for 2 storage box inserts	01185	1	69,99	7
Storage box	for 3 storage box inserts	01142	1	79,99	7
Storage box insert – for small accessories		01148	1	64,99	7
Storage box insert – for microreaction vessels	32 x 0,5 ml	01104	1	64,99	8
Storage box insert – for microreaction vessels	20 x 2 ml	01130	1	64,99	8
Storage box insert – for sample containers and test tubes	20 x Ø 3 - 10 mm	01161	1	64,99	8
Storage box insert – for sample containers and test tubes	15 x Ø 8 - 12 mm	01181	1	64,99	8
Storage box insert – for sample containers and test tubes	12 x Ø 10 - 20 mm	01160	1	64,99	8
Storage box insert – for sample containers and test tubes	8 x Ø 18 - 22 mm	01182	1	64,99	8
Storage box insert – for sample containers and test tubes	6 x Ø 20 - 30 mm	01159	1	64,99	8
Storage box insert – for centrifuge tubes	12 x 15 ml	01128	1	64,99	9
Storage box insert – for centrifuge tubes	6 x 15 ml, 4 x 50 ml	01102	1	64,99	9
Storage box insert – for centrifuge tubes	6 x 50 ml	01129	1	64,99	9
Adapter for LAUDA lid holder		01223	1	49,99	9
AKU® PP-Drinking cup, 500 ml		09006	1	5,99	11
Stainless steel round for insertion into pipette tip boxes		01167	1	2,99	22
Holder for 2 pipettes		01180	1	64,99	10
Holder for 3 Eppendorf pipette holders		01204	1	79,99	10
Holder for 3 pipettes		01145	1	69,99	10
Holder for barcode scanner	suitable for scanners such as Newland NLS-HR22 and Aplic 304674	01163	1	74,99	11
Holder for cup	for AKU® PP-Drinking cup, 500 ml	01213	1	59,99	11
Holder for display	for VESA 75, ball joint	01221	1	229,99	12
Holder for display	for VESA 100, ball joint	01222	1	229,99	12
Holder for display	for VESA 75, fixed	01230	1	169,99	12
Holder for display	for VESA 100, fixed	01231	1	169,99	12
Holder for chromatography column	Ø 5 - 10 mm	01177	1	69,99	12
Holder for chromatography column	Ø 10 - 15 mm	01176	1	69,99	12
Holder for chromatography column	Ø 15 - 25 mm	01165	1	69,99	12
Holder for chromatography column	Ø 25 - 35 mm	01107	1	69,99	12
Holder for chromatography column	Ø 35 - 46 mm	01139	1	69,99	12
Holder for chromatography column	Ø 46 - 56 mm	01179	1	79,99	12

Item	Type	Item no.	PU	Price in EUR	Page
Holder for various KNAUER accessories		01210	1	59,99	13
Holder for Erlenmeyer flasks	Erlenmeyer flasks wide neck 25, 100 ml, Erlenmeyer flasks narrow neck 250, 500 ml	01108	1	69,99	13
Holder for Erlenmeyer flasks	Erlenmeyer flasks wide neck 250 ml, 500 ml, 1000 ml, Erlenmeyer flasks narrow neck 2000 ml, 3000 ml, 5000 ml	01140	1	69,99	13
Holder for glassware with standard taper	ST 14	01111	1	49,99	14
Holder for glassware with standard taper	ST 29	01114	1	49,99	14
Holder for glassware with standard taper plus extension	ST 14	01113	1	59,99	14
Holder for glassware with standard taper plus extension	ST 29	01116	1	59,99	14
Holder for glassware with standard taper plus closure	ST 14	01112	1	64,99	15
Holder for glassware with standard taper plus closure	ST 29	01115	1	64,99	15
Holder for glove / tissue boxes		01143	1	89,99	15
Holder for Inlabtex Serial Dilution Bag box		01216	1	79,99	16
Holder for cannula disposal container	for Vernacare sharpsafe® 11	01174	1	49,99	16
Holder for KNAUER syringe filter box		01212	1	79,99	17
Holder for KNAUER Valve Drive		01217	1	79,99	17
Holder for kitchen roll		01162	1	89,99	18
Holder for cuvettes	10 x (10 mm x 10 mm)	01183	1	64,99	18
Holder for laboratory screw top bottles	10 - 25 ml (GL25)	01120	1	64,99	19
Holder for laboratory screw top bottles	50 ml (GL32)	01121	1	64,99	19
Holder for laboratory screw top bottles	100 - 1000 ml (GL45)	01122	1	64,99	19
Holder for microliter syringes		01155	1	89,99	19
Holder for bin liner	without cover	01202	1	59,99	20
Holder for NETGEAR® ProSAFE GS108		01211	1	69,99	20
Holder for notebook	for 15 inch	01224	1	289,99	21
Holder for pipettes	for Eppendorf Multipipette® E3 / E3x / stram / Xstream	01189	1	59,99	21
Holder for pipettes	for Eppendorf Reference®	01164	1	49,99	21
Holder for pipettes	for Eppendorf Reference® 2	01168	1	49,99	21
Holder for pipettes	for Gilson® PIPETMAN® Classic / G	01166	1	49,99	21
Holder for pipettes	for INTEGRA VIAFLOW / INTEGRA VOYAGER	01184	1	49,99	21
Holder for pipette tip boxes		01147	1	54,99	22
Holder for pipetting aid	for Eppendorf Easytip® 3	01203	1	59,99	22
Holder for pipetting aid	for HIRSCHMANN® pipet- tus®	01188	1	59,99	22

Item	Type	Item no.	PU	Price in EUR	Page
Holder for pipetting aid	for INTEGRA-Pipetboy 2	01154	1	59,99	22
Holder for pipetting aid	for INTEGRA-Pipetboy pro	01156	1	59,99	22
Holder for pipetting aid	for BRAND accu-jet® pro	01178	1	59,99	22
Holder for pipetting aid	for BRAND accu-jet® S	01157	1	59,99	22
Holder for pipetting aid	for Thermo Scientific Matrix	01186	1	59,99	22
Holder for pipetting aid	for TPP Turbo-Fix	01207	1	59,99	22
Holder for pipetting aid	for VWR® Powerpette Pro	01199	1	59,99	22
Holder for separating funnel / shaking funnel	Ø 45 - 85 mm	01152	1	59,99	23
Holder for separating funnel / shaking funnel	Ø 85 - 105 mm	01149	1	59,99	23
Holder for separating funnel / shaking funnel	Ø > 105 mm	01153	1	59,99	23
Holder for sensor systems	for essentim sphere	01208	1	69,99	23
Holder for sensor systems	for essentim spot	01209	1	69,99	23
Holder for sensor systems	for essentim spider	01205	1	79,99	23
Holder for open fronted storage bins	suitable for storage bins such as SCHÄFER LF110, LF210, LF211, LF220	01158	1	69,99	24
Holder for wash bottles	for 500 ml (up to ø 75 mm)	01146	1	59,99	24
Holder for wash bottles	for 1000 ml (up to ø 95 mm)	01173	1	69,99	24
Holder for tablet		01124	1	99,99	25
Holder for well plates		01125	1	54,99	25
Holder for ZEBRA printer ZD621		01201	1	179,99	26
Holder for ZEBRA hand-held scanner TC22 standing at an angle		01206	1	79,99	26
Holder for draining	6x drainage prong, I-shaped	01225	1	59,99	27
Holder for draining	6x drainage prong, Y-shaped	01226	1	69,99	27
Holder for draining	3x drainage prong, I-shaped and 3x drainage prong, Y-shaped	01227	1	64,99	27
Shelf	Width 299 mm, Depth 185 mm	01193	1	119,99	27
Dispenser for adhesive tapes		01172	1	149,99	28
Dispenser for Parafilm® rolls	without tear-off blade	01171	1	149,99	28
Dispenser for Parafilm® rolls	with tear-off blade	01214	1	199,99	28
Dispenser for glass pasteur pipettes	150 mm	01169	1	149,99	29
Dispenser for glass pasteur pipettes	230 mm	01170	1	169,99	29
Stand rod adapter	ø 10 mm	01218	1	79,99	29
Stand rod adapter	ø 12 mm	01126	1	79,99	29
Additional core for adhesive tape dispenser	approx. 25 mm (1 inch) / approx. 51 mm (2 inch)	01198	1	14,99	28
Additional core for adhesive tape dispenser	approx. 25 mm (1 inch) / approx. 76 mm (3 inch)	01175	1	24,99	28

Additional Accessories

Item	Type	Item no.	PU	Price in EUR	Page
Lid for storage box	for 3 storage box inserts	01229	1	34,99	30
Lid for storage box	for 2 storage box inserts	01228	1	49,99	30
Lid for bin liner holder		01192	1	19,99	30

SmartStand

Item	Type	Item no.	PU	Price in EUR	Page
SmartStand for well plates		09001	1	129,99	31

Technical Information

Technical information:

- All products are washable with water and mild detergent.
- In addition all components can be wiped with alcohol-based disinfectants.
- The plastics and metals used are 100% recyclable.
- Our products are autoclavable up to 134 °C and 3 bar.
- The plastic elements are resistant to heat deformation under normal laboratory conditions.
- To gently clean the anodised aluminium rail, we recommend using a pH-neutral hand-washing detergent. Ensure that the hand wash detergent is odourless and contains no perfumes, greases or additives. Clean the slot rail by hand and rinse thoroughly with clean water.

Note on joint tolerances for ST14:

- Maximum outside diameter at the narrowest point of the cone joint: 17.2 mm (larger diameters are too wide for the holder and cannot be held)
- Minimum outer diameter at the widest point of the cone joint: 18.0 mm (smaller diameters are too narrow for the holder and would fall out)

Note on joint tolerances for ST29:

- Maximum outside diameter at the narrowest point of the cone joint: 31.4 mm (larger diameters are too wide for the holder and cannot be held)
- Minimum outer diameter at the widest point of the cone joint: 32.8 mm (smaller diameters are too narrow for the holder and would fall out)

IMPORTANT TO NOTE: The products of the SmartIntegrate line, the SmartRack® and SmartRack® mini as well as the SmartCase and SmartStand are not suitable for the permanent storage of hazardous substances.

We reserve the right to make technical changes, material or model changes, errors and misprints excepted. If you have any further questions, please do not hesitate to contact us at contact@better-basics.de.

For more detailed information, please visit our website:

Functionality of SLS 3D printing:	www.sls-3d-printing.better-basics-laborbedarf.de
Resistance to UV light:	www.uv-resistance.better-basics-laborbedarf.de
Resistance to aggressive media:	www.material-resistance.better-basics-laborbedarf.de
Sustainability and design-efficient manufacturing:	www.sustainability.better-basics-laborbedarf.de



Photo: Building of the Better Basics Laborbedarf GmbH in Dresden

YOU ARE INTERESTED IN OUR PRODUCTS OR YOU HAVE FURTHER QUESTIONS?

Get in touch with us:

contact@better-basics.de

You can also find more information and news about our product innovations at

www.Better-Basics-Laborbedarf.de

Better Basics Laborbedarf GmbH

Loebtauer Str. 69

01159 Dresden

Germany

T: +49 (0) 178 9617 577 (Direct customers)

T: +49 (0) 176 6233 8026 (Sales partners)

E: contact@better-basics.de

W: www.Better-Basics-Laborbedarf.de

Our request to you

Our products such as our SmartRack® and our SmartIntegrate products are only as good as their users. With our roots in the Leibniz Institute for Polymer Research Dresden, our company Better Basics Laborbedarf stands for the tradition of “Made in Germany” quality. We feel in multiple aspects committed to this quality label for German companies, which stands for high innovative strength. Even beyond our product catalogue, the same applies for us at Better Basics: Just contact us if you have ideas for extensions or improvements or if you notice something in your daily laboratory routine that we should optimise right away.

We want to create the best possible working environment for you in the laboratory and we are grateful for any comments and inspiration. We can jointly shape the future of laboratory work by cooperating with scientists and researchers. Just like you, we’re always looking for a way to make our world a little better every day.

Impress:

Better Basics Laborbedarf GmbH
Loebtauer Str. 69, G-01159 Dresden

E-Mail: info@better-basics.de

CEO: Mario Schneider

Commercial register: HRB 39524; Register court: Amtsgericht Dresden

VAT registration number: DE327508812

