

SmartIntegrate

Integrative organisational solutions from
Better Basics for Wesemann laboratory equipment

SmartFrame - Fume hood

SmartEdge - Fume hood

SmartRail - Workbench

Contact us for a personalised consultation:
Contact@Better-Basics.de



SmartFrame elite WE 001 in a
Wesemann Delta 30 fume hood
Click [HERE](#) for the product video.

The SmartRail WE 001 and SmartRail WE 002 have been specifically developed for use on Wesemann media channels with functional slots. As an extension of the SmartRail WE 001, the SmartRail WE 002 includes an additional rotary lock to enable stable mounting on the front functional slot of the media channels. The construction consists of a black anodised aluminium slot rail with side brackets made of PA12 and stainless steel, featuring a practical snap-in function to ensure secure attachment to the front functional slot of the media channels.

The versatile SmartRail allows flexible placement options within the laboratory. A larger number of SmartRails provides multiple possibilities for module placement and enables the individual adaptation of the laboratory environment to specific requirements. All Better Basics modules are fully compatible with the SmartRail.

Click [HERE](#) for the product video.

SmartRail

Solutions for the laboratory workbench



SmartRail WE 002 in a
Wesemann Delta 30 workbench



SmartEdge elite WE 001 in a
Wesemann Delta 30 fume hood

SmartEdge

Solutions for the fume hood

The SmartEdge elite WE 001 for Wesemann Delta 30 fume hoods consists of two stainless-steel L-profiles, which are mounted on the left and right onto the stand receptacles inside the fume hood. Depending on the width of the fume hood, a variable number of slot rails can be flexibly inserted into these L-profiles. This creates additional storage space and a stable, adaptable organisational base that ensures laboratory materials are safely, conveniently, and clearly arranged within the fume hood. A larger number of slot rails provides multiple options for module placement and enables the individual adaptation of the laboratory environment to specific requirements. All Better Basics modules are fully compatible with the SmartEdge.

Click [HERE](#) for the product video.

The SmartFrame elite WE 001 and WE 002 have been developed for use in the Wesemann Delta 30 laboratory fume hood. They consist of a brushed stainless-steel frame with integrated slot rails.

Every third slot rail mount on both sides of the SmartFrame elite features a recessed front section, which facilitates the precise horizontal alignment of the slot rails. All Better Basics modules are fully compatible with the SmartFrame.

Click [HERE](#) for the product video.

SmartFrame

Solutions for the fume hood



SmartFrame elite WE 002 in a
Wesemann Delta 30 fume hood

Over 100 modules

One system for endless possibilities to optimise every workplace and workflow



All Better Basics modules that can be directly mounted onto slot rails feature an integrated click mechanism, providing users with tactile feedback when positioning the holding modules within the slot rail. Since all holding modules can be freely positioned along a slot rail, you have complete flexibility in arranging them to suit your specific workflow. Once clicked into place, the modules remain securely positioned, even when the slot rail itself is moved to another working level. This ensures that both you and your work are protected against sample loss and workplace accidents.

Item numbers and price

Your guide to ordering and calculation

As of October 2025

Item	Type	PU	Item-N°	Price in EUR
SmartEdge elite WE 001	elite	1	04011	229,99
Slot rail for SmartEdge WE 001	765 mm	1	04334	59,99
Slot rail for SmartEdge WE 001	865 mm	1	04338	59,99
Slot rail for SmartEdge WE 001	915 mm	1	04336	64,99
Slot rail for SmartEdge WE 001	1065 mm	1	04335	69,99
Slot rail for SmartEdge WE 001	1265 mm	1	04339	74,99
Slot rail for SmartEdge WE 001	1365 mm	1	04337	79,99
Slot rail for SmartEdge WE 001	1665 mm	1	04340	89,99
SmartRail 300 WE 002	300 mm	1	03038	79,99
SmartRail 600 WE 002	600 mm	1	03039	99,99
SmartRail 900 WE 002	900 mm	1	03040	119,99
SmartRail 1200 WE 002	1200 mm	1	03041	139,99
SmartRail 1500 WE 002	1500 mm	1	03042	159,99
SmartFrame elite WE 001	elite	1	05020	529,00
SmartFrame elite WE 002	elite	1	05021	579,00

Better Basics Laborbedarf - Unlock Your Lab Potential



Are you interested in our products or do you have further questions?

Get in touch with us: Contact@Better-Basics.de

You can also find more information and news about our product innovations at: www.Better-Basics-Laborbedarf.de

Better Basics Laborbedarf GmbH | Loebtauer Str. 69
01159 Dresden | Germany

T: +49 (0) 176 6233 8026 (Sales partner)

T: +49 (0) 178 9617 577 (direct customers)

Better Basics Laborbedarf

is a leading provider of innovative organisational systems such as the SmartRack® and the SmartIntegrate product line, revolutionising laboratory workstations. Our award-winning solutions optimise efficiency, safety, and ergonomics. By improving space utilisation and offering modular adaptability, they reduce walking distances and increase satisfaction in the laboratory. As a reliable partner, we develop customised systems tailored to individual requirements – flexible and adaptable, wherever you need them.

Impressum:

Better Basics Laborbedarf GmbH | Löbtauer Str. 69 | G-01159 Dresden
E-Mail: Info@Better-Basics.de | Geschäftsführer: Mario Schneider | Handelsregister: HRB 39524
Registergericht: Amtsgericht Dresden | Ust.-Id.-Nr.: DE327508812